

FLEXIBLE SOLUTIONS FOR TRANSPORT



EQUIPMENT LOADING RAMPS

CONVENIENT EQUIPMENT LOADING RAMPS ARE MADE OF STURDY LIGHTWEIGHT ALUMINUM. THEY ENABLE FAST AND EASY LOADING OF THE PORTABLE X-RAY MACHINE INTO VANS AND CARS. THESE RAMPS FOLD IN HALF (36"), AND HAVE CARRY HANDLES WELDED ON THEM FOR EASY CARRYING AND STORAGE. STRONG CONSTRUCTION WITH A 600LB WEIGHT LOAD CAPACITY.

AN OPTIONAL SPACER BAR IS AVAILABLE ENSURING INCREASED STABILITY FOR VAN ACCESS OR PARTICULARLY STEEP APPLICATIONS. THIS BAR CAN BE MOUNTED ON THE THRESHOLD OF A VAN DOORWAY, AND WILL SECURELY LOCK EACH TRACK IN PLACE.



ORDERING INFORMATION

P/N 420059: LOADING RAMPS, 6' FOLDING SET

P/N 420060: BAR, RAMP SPACER (OPTIONAL)

CONVENIENT FOLDING HINGE & CARRY HANDLES

