

# Mobile DR Chest **Holder**

## Flexible Solutions for Urgent Care



### Key Advantages

- ◆ Radiographic Area: 14" x 17", or 17 x 14"
- ◆ Panel can be positioned for portrait or landscape
- ◆ Vertical travel: 6" to 74"
- ◆ Pushbutton counterbalanced pulley system
- ◆ Removable grid: 103 line, 6:1, 40"-72" SID
- ◆ User friendly
- ◆ Hospital grade casters
- ◆ Maximum Dimensions
  - Height: 78.5"
  - Width: 24.0"
  - Depth: 22.0"
  - Weight: 87 lbs

