

Mobile X-Ray Shield

Flexible Solutions for Urgent Care



- ◆ Reliable radiation protection for radiology (30" x 72").
- ◆ Ribbed construction for strength and rigidity.
- ◆ 1.55 mm lead liner centered between two steel walls for a total thickness of 7/8"
- ◆ Built in 8" x 10" lead glass window
 - ◆ (1.7 mm minimum, lead equivalency).
- ◆ Window centered in the panel with its center located 60" above the floor.
- ◆ Shield extends to 0.25" above floor.
- ◆ Stability and ease of movement are assured by the 3/16" thick steel legs, with ball bearing casters.
- ◆ Minor assembly required
- ◆ Ships in protective wood crate.
- ◆ Made in USA