

Table Top DR Panel Holder



Taking shoot-through exams was difficult enough with conventional radiography, now with the advent of mobile DR the image receptor has become larger, more cumbersome and increasingly heavy.

The DRP Table Top Holder for Varian's 4336R allows the radiographer to position the mobile DR panel with ease and precision.

The DRP-TTH can also be angled caudal and cephalic depending on the examination.

Product Advantages

- ◆ Rotates 360 degrees for caudal and cephalic angling
- ◆ Extends 27" high
- ◆ Size adjustable
- ◆ Long lasting parts made in the USA